




Plant actin Monoclonal Antibody

| | |
|----------------------------|--|
| Product Code | CSB-MA000307 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Recombinant Protein |
| Raised In | Mouse |
| Species Reactivity | Various Plants |
| Tested Applications | ELISA, WB; Recommended dilution: WB: 1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Various Plants |

Image



Western blot analysis of Arabidopsis, diluted at
1) 1:5000 2) 1:10000