




TF Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA000313 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02787 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | TF; Serotransferrin; Transferrin; Beta-1 metal-binding globulin; Siderophilin |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Gene Names | TF |

Image



Western blot analysis of Human serum, mAb diluted at 1:2000.