




FH Monoclonal Antibody

Product Code	CSB-MA000317
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07954
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,IF;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Fumarate hydratase, mitochondrial
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	FH
Clone No.	20C1

Image



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.