




FN1 Monoclonal Antibody

Product Code	CSB-MA000318
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P02751
Immunogen	Synthetic Peptide
Raised In	Mouse
Species Reactivity	Mouse,Rat,Human
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	FN1; FN; Fibronectin; FN; Cold-insoluble globulin; CIG
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Gene Names	FN1

Image



Western blot analysis of 1) Mouse Brain Tissue, 2) Rat Brain tissue, diluted at 1:2000.