



LTF Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA000322 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P02788 |
| Immunogen | Synthetic Peptide |
| Raised In | Mouse |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |
| Alias | LTF; LF; Lactotransferrin; Lactoferrin; Talalactoferrin |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | LTF |
| Clone No. | 5C5 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |