




IgG (Fc specific, Heavy Chain) Monoclonal Antibody

Product Code	CSB-MA000329
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	Rabbit
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Monoclonal Antibody
Immunogen Species	Oryctolagus cuniculus (Rabbit)

Image



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.