




# Flag-Tag Antibody

<b>Product Code</b>	CSB-PA000337
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Tested Applications</b>	WB;WB:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG

## Image



1µg Flag fusion protein+ Primary antibody dilution at 1) 1:5000 2) 1:10000 Transfected 293 cell lysate with Flag-tag fusion protein+ Primary antibody dilution at 1:10000