




# HA-Tag Antibody

<b>Product Code</b>	CSB-PA000338
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Tested Applications</b>	WB;WB:1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Isotype</b>	IgG
<b>Alias</b>	0

**Image**

2 $\mu$ g HA fusion protein+ Primary antibody dilution at 1) 1:5000 2) 1:10000

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.