




GST-Tag Antibody

Product Code	CSB-PA000341
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Tested Applications	WB;WB:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG

Image



0.5µg GST fusion protein+ Primary antibody dilution at 1) 1:5000 2) 1:10000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.