





ACTB Antibody

Product Code	CSB-PA000350
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat, Rabbit, Chicken, Monkey, Sheep
Tested Applications	WB, IHC; WB: 1:3000 IHC: 1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB

Image



IHC staining of Human ovary tissue, diluted at 1:200.



Western blot analysis of HepG2 (1), Rat liver (2), Mouse kidney (3), Rabbit testis (4), Sheep lung (5), 293T (6), diluted at 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.