





ACTB Antibody

Product Code	CSB-PA000350
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P60709
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat,Rabbit,Chicken,Monkey,Sheep,X
Tested Applications	WB,IHC;WB:1:3000IHC:1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	ACTB; Actin, cytoplasmic 1; Beta-actin
Immunogen Species	Homo sapiens (Human)
Target Names	ACTB

Image



IHC staining of Human ovary tissue, diluted at 1:200.



Western blot analysis of HepG2 (1), Rat liver (2), Mouse kidney (3), Rabbit testis (4), Sheep lung (5), 293T (6), diluted at 1:5000.