




Plant-actin Antibody

Product Code	CSB-PA000352
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Arabidopsis Nicotiana tabacum Oryza sativa Hordeum vulgare
Tested Applications	WB;WB:1:3000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0
Immunogen Species	Arabidopsis Nicotiana tabacum Oryza sativa Hordeum vulgare

Image



Western blot analysis of Nicotiana tabacum leaf (1), Nicotiana tabacum root (2), Nicotiana seed (3), Oryza sativa flower (4), Oryza sativa seed (5), Oryza sativa leaf (6), Hordeum vulgare seed (7), Hordeum vulgare leaves (8), diluted at 1:3000.