




ECFP-Tag Polyclonal Antibody

Product Code	CSB-PA000355
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Mouse
Species Reactivity	N/A
Tested Applications	WB; WB: 1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0

Image

Western blot analysis of ECFP Recombinant protein, diluted at 1) 1:5000 2) 1:10000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.