




Plant GST Antibody

Product Code	CSB-PA000356
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Nicotiana tabacum
Tested Applications	WB;WB:1:3000-10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Immunogen Species	Nicotiana tabacum (Common tobacco)

Image



10µg total Nicotiana tabacum lysate + primary antibody dilute at 1:3000

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.