




EFHD1 Antibody

Product Code	CSB-PA000358
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BUP0
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB;WB:1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	EF-hand domain-containing protein D1
Immunogen Species	Homo sapiens (Human)
Target Names	EFHD1

Image



Western blot analysis of 293T cell lysate, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.