




FH Antibody

Product Code	CSB-PA000359
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07954
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	Fumarate hydratase, mitochondrial
Immunogen Species	Homo sapiens (Human)
Target Names	FH

Image



Western blot analysis of 1) Hela, 2) Mouse Brain, 3) HepG2, diluted at 1:2000.