




PTGS1 Antibody

Product Code	CSB-PA000360
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23219
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB;WB:1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	PTGS1; COX1; Prostaglandin G/H synthase 1; Cyclooxygenase-1; COX-1; Prostaglandin H2 synthase 1; PGH synthase 1; PGHS-1; PHS 1; Prostaglandin-endoperoxide synthase 1
Immunogen Species	Homo sapiens (Human)
Target Names	PTGS1

Image



Western blot analysis of HeLa, diluted at 1:1000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.