




# PTGS2 Antibody

<b>Product Code</b>	CSB-PA000361
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P35354
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	WB;WB:1:1000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
<b>Isotype</b>	IgG
<b>Alias</b>	PTGS2; COX2; Prostaglandin G/H synthase 2; Cyclooxygenase-2; COX-2; PHS II; Prostaglandin H2 synthase 2; PGH synthase 2; PGHS-2; Prostaglandin-endoperoxide synthase 2
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	PTGS2

## Image



Western blot analysis of HeLa, diluted at 1:1000.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.