




Alliinase Antibody

Product Code	CSB-PA000362
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified Protein
Raised In	Rabbit
Species Reactivity	Allium sativum L
Tested Applications	WB;WB:1:3000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0
Immunogen Species	Allium sativum L

Image



Western blot analysis of purified alliinase, diluted at 1:5000