




GAPDH Antibody

Product Code	CSB-PA000364
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04406
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	WB;WB:1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH
Immunogen Species	Homo sapiens (Human)
Target Names	GAPDH

Image



Western blot analysis of Zebrafish skeletal muscle, diluted at 1:5000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.