




mOrange Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA000367 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Recombinant Protein |
| Raised In | Rabbit |
| Tested Applications | WB;WB:1:2000-5000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol. |
| Isotype | IgG |

Image

Western blot analysis of mOrange recombinant protein, diluted at 1:5000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.