



BFP Polyclonal Antibody

Product Code	CSB-PA000369
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	N/A
Tested Applications	WB; WB: 1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide
Isotype	IgG
Product Type	Polyclonal Antibody
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.