




Rubisco Polyclonal Antibody

Product Code	CSB-PA000372
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified Protein
Raised In	Mouse
Species Reactivity	Various Plants
Tested Applications	WB; WB: 1:3000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	0

Image



Western blot analysis of Arabidopsis thaliana, diluted at 1:5000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.