




BFP Antibody

Product Code	CSB-PA000375
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Tested Applications	WB;WB:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Polyclonal Antibody

Image



Western blot analysis of EBFP recombinant protein, diluted at 1:5000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.