




CREB1 Antibody

Product Code	CSB-PA000376
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16220
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	CREB1; Cyclic AMP-responsive element-binding protein 1; CREB-1; cAMP-responsive element-binding protein 1
Immunogen Species	Homo sapiens (Human)
Target Names	CREB1

Image



Western blot analysis of HeLa, Rabbit pAb diluted at 1:2000.