




LMNB1 Antibody

Product Code	CSB-PA000377
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20700
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000-3000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	LMNB1; LMN2; LMNB; Lamin-B1
Immunogen Species	Homo sapiens (Human)
Target Names	LMNB1

Image



Western blot analysis of 1) 293T, 2) Mouse Brain, 3) Rat Brain tissue, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.