




VCL Antibody

Product Code	CSB-PA000379
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18206
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	VCL; Vinculin; Metavinculin
Immunogen Species	Homo sapiens (Human)
Target Names	VCL

Image



Western blot analysis of 1) Hela, 2) K562, 3) Mouse Brain, 4) Rat Brain tissue, diluted at 1:1000.