




Histone H1 Antibody

Product Code	CSB-PA000381
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02539/P16403/P16402/P10412
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB; WB: 1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4; Histone H1b; Histone H1s-4
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E

Image



Western blot analysis of 1) HeLa, 2) 3T3, diluted at 1:1000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.