




ACTA1 Antibody

Product Code	CSB-PA000383
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68133
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	ACTA1; ACTA; Actin, alpha skeletal muscle; Alpha-actin-1
Immunogen Species	Homo sapiens (Human)
Target Names	ACTA1

Image



Western blot analysis of 1) HeLa, 2) 3T3, 3) Rat Kidney tissue, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.