




# GAP43 Antibody

<b>Product Code</b>	CSB-PA000387
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P17677
<b>Immunogen</b>	Recombinant Protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB;WB:1:1000-3000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	GAP43; Neuromodulin; Axonal membrane protein GAP-43; Growth-associated protein 43; Neural phosphoprotein B-50; pp46
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	GAP43

## Image



Western blot analysis of 1) HeLa, 2) Mouse Brain, 3) Rat Brain tissue, diluted at 1:3000.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.