




MYC Antibody

Product Code	CSB-PA000389
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01106
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB;WB:1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Alias	MYC; BHLHE39; Myc proto-oncogene protein; Class E basic helix-loop-helix protein 39; bHLHe39; Proto-oncogene c-Myc; Transcription factor p64
Immunogen Species	Homo sapiens (Human)
Target Names	MYC

Image



Western blot analysis of HeLa, Rabbit pAb diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.