




BSA Antibody

Product Code	CSB-PA000390
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Cow
Tested Applications	WB,ELISA;WB:1:2000-5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Cow

Image



Western blot analysis of Bovine Serum Albumin, Rabbit pAb diluted at 1) 1:2000 2) 1:5000

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.