




KRT18 Antibody

Product Code	CSB-PA000392
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05783
Immunogen	Recombinant Protein
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	KRT18; CYK18; PIG46; Keratin, type I cytoskeletal 18; Cell proliferation-inducing gene 46 protein; Cytokeratin-18; CK-18; Keratin-18; K18
Immunogen Species	Homo sapiens (Human)
Target Names	KRT18

Image



Western blot analysis of HeLa, Rabbit pAb diluted at 1:2000.

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.