




Acetyl-Histone H3 (Lys27) Antibody

Product Code	CSB-PA000410
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P68431/Q71DI3/P84243
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:200-500
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)

Image



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.