




Mono-Methyl-Histone H2B (Lys5) Antibody

Product Code	CSB-PA000426
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A08/P33778/P62807
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H2BI; H2BFK
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H2BA/HIST1H2BB/HIST1H2BC

Image



Western blot analysis of 1) HeLa, 2) 3T3, diluted at 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.