




# Acetyl-Histone H2A (Lys5) Antibody

<b>Product Code</b>	CSB-PA000435
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P0C0S8/Q6F113/Q7L7L0
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB;WB:1:1000-2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	HIST1H2AG; H2AFP; HIST1H2AI; H2AFC; HIST1H2AK; H2AFD; HIST1H2AL; H2AFI; HIST1H2AM; H2AFN; Histone H2A type 1; H2A.1; Histone H2A/p; HIST2H2AA3; H2AFO; HIST2H2AA; HIST2H2AA4; Histone H2A type 2-A; Histone H2A.2; Histone H2A/o; HIST3H2A; Histone H2A type 3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	HIST1H2AG/HIST2H2AA3/HIST3H2A

## Image



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.