




Acetyl-Histone H2A.Z (Lys7) Antibody

Product Code	CSB-PA000438
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0S5
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:1000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	H2AFZ; H2AZ; Histone H2A.Z; H2A/z
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFZ

Image



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.