




Phospho-Histone H1 (Ser1) Antibody

Product Code	CSB-PA000447
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02539/P16403/P16402/P10412
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	HIST1H1A; H1F1; Histone H1.1; Histone H1a; HIST1H1C; H1F2; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1; HIST1H1D; H1F3; Histone H1.3; Histone H1c; Histone H1s-2; HIST1H1E; H1F4; Histone H1.4; Histone H1b; Histone H1s-4
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H1A/HIST1H1C/HIST1H1D/HIST1H1E

Image



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.