




# Phospho-Histone H2A.X (Ser139) Antibody

<b>Product Code</b>	CSB-PA000451
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P16104
<b>Immunogen</b>	Synthetic Peptide
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human,Mouse,Rat
<b>Tested Applications</b>	WB;WB:1:1000-2000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Isotype</b>	IgG
<b>Alias</b>	H2AFX; H2AX; Histone H2A.x; H2a/x
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	H2AFX

## Image



Western blot analysis of extracts from HeLa cells, 1:2000.