




Acetyl-Histone H2A.X (Lys5) Antibody

Product Code	CSB-PA000455
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P16104
Immunogen	Synthetic Peptide
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Alias	H2AFX; H2AX; Histone H2A.x; H2a/x
Immunogen Species	Homo sapiens (Human)
Target Names	H2AFX

Image



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.