




Acetyl-Histone H2B (Lys5) Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA000457 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96A08/P33778/P62807 |
| Immunogen | Synthetic Peptide |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB;WB:1:1000-2000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Isotype | IgG |
| Alias | HIST1H2BA; TSH2B; Histone H2B type 1-A; Histone H2B, testis; Testis-specific histone H2B; HIST1H2BB; H2BFF; Histone H2B type 1-B; Histone H2B.1; Histone H2B.f; H2B/f; HIST1H2BC; H2BFL; HIST1H2BE; H2BFH; HIST1H2BF; H2BFG; HIST1H2BG; H2BFA; HIST1H2BI; H2BFK |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | HIST1H2BA/HIST1H2BB/HIST1H2BC |

Image



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:5000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.