




Mono-Methyl-Histone H4 (Arg3) Antibody

Product Code	CSB-PA000463
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62805
Immunogen	Synthetic Peptide
Raised In	mouse
Species Reactivity	Human
Tested Applications	WB;WB:1:1000-2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Target Names	HIST1H4A

Image



Western blot analysis of extracts from HeLa cells, 1:2000.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.