




Phospho-BAD (S155) Antibody

Product Code	CSB-PA000483
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92934
Immunogen	Synthesized peptide derived from Human Bad around the phosphorylation site of S155.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BAD; BBC6; BCL2L8; Bcl2 antagonist of cell death; BAD; Bcl-2-binding component 6; Bcl-2-like protein 8; Bcl2-L-8; Bcl-XL/Bcl-2-associated death promoter
Immunogen Species	Homo sapiens (Human)
Target Names	BAD

Image



Western Blot analysis of 293 cells using Phospho-Bad (S155) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.