




Phospho-BCR (Y177) Antibody

Product Code	CSB-PA000486
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11274
Immunogen	Synthesized peptide derived from Human Bcr around the phosphorylation site of Y177.
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	BCR; BCR1; D22S11; Breakpoint cluster region protein; Renal carcinoma antigen NY-REN-26
Immunogen Species	Homo sapiens (Human)
Target Names	BCR

Image



Western Blot analysis of HuvEc cells using Phospho-Bcr (Y177) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.