




Phospho-CEBPB (T235) Antibody

Product Code	CSB-PA000487
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17676
Immunogen	Synthesized peptide derived from Human C/EBP β around the phosphorylation site of T235.
Raised In	Rabbit
Species Reactivity	Human, Mouse, Rat
Tested Applications	WB, IHC, IF, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CEBPB; LAP; TCF5; PP9092; CCAAT/enhancer-binding protein beta; C/EBP beta; Liver activator protein; Nuclear factor NF-IL6; Transcription factor 5; TCF-5
Immunogen Species	Homo sapiens (Human)
Target Names	CEBPB

Image



Western Blot analysis of HepG2 cells using Phospho-C/EBP β (T235) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.