





Phospho-CAV1 (Y14) Antibody

Product Code	CSB-PA000491
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q03135
Immunogen	Synthesized peptide derived from Human Caveolin-1 around the phosphorylation site of Y14.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, ELISA; WB:1:500-1:2000, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CAV1; CAV; Caveolin-1
Immunogen Species	Homo sapiens (Human)
Target Names	CAV1

Image



Western Blot analysis of HeLa cells using Phospho-Caveolin-1 (Y14) Polyclonal Antibody



Western Blot analysis of A549 cells using Phospho-Caveolin-1 (Y14) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.