




Phospho-KIT (Y721) Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA000502 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P10721 |
| Immunogen | Synthesized peptide derived from Human c-Kit around the phosphorylation site of Y721. |
| Raised In | Rabbit |
| Species Reactivity | Human,Mouse,Rat |
| Tested Applications | WB, ELISA; WB:1:500-1:2000, ELISA:1:10000 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Purification Method | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Isotype | IgG |
| Alias | KIT; SCFR; Mast/stem cell growth factor receptor Kit; SCFR; Piebald trait protein; PBT; Proto-oncogene c-Kit; Tyrosine-protein kinase Kit; p145 c-kit; v-kit Hardy-Zuckerman 4 feline sarcoma viral oncogene homolog; CD antigen CD117 |
| Immunogen Species | Homo sapiens (Human) |
| Target Names | KIT |

Image



Western Blot analysis of HepG2 cells using Phospho-c-Kit (Y721) Polyclonal Antibody

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.