





Phospho-GJA1 (S368) Antibody

Product Code	CSB-PA000504
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17302
Immunogen	Synthesized peptide derived from Human Connexin 43 around the phosphorylation site of S368.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:20000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	GJA1; GJAL; Gap junction alpha-1 protein; Connexin-43; Cx43; Gap junction 43 kDa heart protein
Immunogen Species	Homo sapiens (Human)
Target Names	GJA1

Image



Western Blot analysis of HuvEc cells using Phospho-Connexin 43 (S368) Polyclonal Antibody



Western Blot analysis of A549 cells using Phospho-Connexin 43 (S368) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.