




Phospho-CCNB1 (S126) Antibody

Product Code	CSB-PA000507
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14635
Immunogen	Synthesized peptide derived from Human Cyclin B1 around the phosphorylation site of S126.
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCNB1; CCNB; G2/mitotic-specific cyclin-B1
Immunogen Species	Homo sapiens (Human)
Target Names	CCNB1

Image



Western Blot analysis of HepG2 cells using Phospho-Cyclin B1 (S126) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.