





Phospho-CCNE1 (T395) Antibody

Product Code	CSB-PA000508
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24864
Immunogen	Synthesized peptide derived from Human Cyclin E1 around the phosphorylation site of T395.
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	WB, IHC, ELISA; WB:1:500-1:2000, IHC:1:100-1:300, ELISA:1:10000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Isotype	IgG
Alias	CCNE1; CCNE; G1/S-specific cyclin-E1
Immunogen Species	Homo sapiens (Human)
Target Names	CCNE1

Image



Western Blot analysis of HeLa cells using Phospho-Cyclin E1 (T395) Polyclonal Antibody



Western Blot analysis of A549 cells using Phospho-Cyclin E1 (T395) Polyclonal Antibody

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.